

No One Knows What Life Is

The greatest scientists believe it
may be a form of electricity.

The results of

Violet Ray Treatment

almost warrant the belief in the
age of miracles. It shines forth
as the outstanding method of
achieving the miraculous.

It is estimated that in this country two hundred and fifty million dollars are lost every year, through sicknesses not grave enough to demand the services of a physician; while the loss that results from over-fatigue or the lowering of vitality and endurance is, from the economic point of view, even greater than that caused by illness.

Many of the ablest specialists in electrical treatment contend that 90 per cent of this loss could be prevented by the proper use of the VIOLET RAY.

Health Vigor and Strength



*How the
Star-Rite Violet Ray
increases Vitality and
drives away Disease*

FITZGERALD MFG. CO., TORRINGTON, CONN., U. S. A.

beneficial, and with treatment of this kind for five to six minutes daily, relief should be felt.

PILES—VIOLET RAY is recommended for the treatment of this affliction. Electrode No. 8, properly lubricated before insertion, is recommended. We suggest that in this treatment you confer with your physician.

RHEUMATISM—For rheumatism, strong currents are recommended. Electrode No. 1 is recommended and should be applied to the affected area until the pain lessens. An application of five to ten minutes is advisable, and can be applied as often as the user desires. We suggest treatment through a cloth, or through layers of thick Turkish toweling.

SCALP TREATMENT—See **ALOPECIA**, **DANDRUFF** and **GRAY HAIR**.

SKIN DISEASES—See **ACNE** and **ECZEMA**.

SCIATICA, **STIFF JOINTS** and **STIFF NECK**—See **RHEUMATISM**.

SORE THROAT—Electrode No. 3 is used with a mild current. Keep in close contact with the skin and use for three to five minutes. If it is a stubborn case, use twice daily. Applications of cloths, wrung out of cold water, may be employed after each treatment, leaving the cloth on for several hours if necessary.

SPRAINS—Contact applications with Electrode No. 1, using a mild current with now and then sharp sparks through a towel, as a rule give considerable relief.

ULCERS—See **ABSCESSES** and **ACNE**.

WARTS—Use Electrode No. 1, letting the sparks pass through the side instead of the face of the electrode to the affected spot. Use as strong a current as possible. Treatment should last three to five seconds. Apply every three days.


GUARANTY

THE STAR-RITE VIOLET RAY is guaranteed for one year from date of purchase, against any defects in material or workmanship of any kind. **THE STAR-RITE VIOLET RAY** is made out of the best material, selected for its efficiency and durability, and the machines manufactured and assembled by skilled workmanship.

TORRINGTON, CONNECTICUT


MODEL A

**\$12.50
Complete**




MODEL B

**\$25.00
Complete**



MODEL C

**\$50.00
Complete**



STAR-RITE VITALATOR

**\$12.50
Complete**

